



Use the F-750 to measure quality indicators in liquids, powders, and irregular solids with the F-750LK Liquids and Powders Kit. The table-mount Liquids Kit enables the operator to measure liquids in the included high-quality glass cuvettes (2 x 1 cm path length, 1 x 5 cm path length) which fit down inside the adjustable holder. For powders, irregular solids, and viscous fluids, use the custom-built petri dish holder to take measurements in one of the 5 included plastic petri dishes.



Excellent for use with a variety of liquids, powders, irregular solids such as ground meats or grains, as well as viscous fluids like peanut butter, molasses, and cooking oils.

